

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	Den2-Den4
Description:	A fusion protein containing dengue subtype 2 and 4 partial envelope sequence is expressed from <i>E. coli</i> , each part in this fusion protein contains 170 amino acids from residue 46 to 217. This protein is fused with 6-His fusion partner.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>95% pure (10% SDS-PAGE Coomassie blue staining). Proprietary chromatographic technique. (Re-check on 10% SDS-PAGE on Oct. 16, 2023, show good pattern.)
Concentration:	1.46mg/ml (Lot 07/2014)
Buffer:	Phosphate buffered saline, pH 7.4
Preservative:	0.05% Sodium azide
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

