


BioassayPlus Inc

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	Den3-N
Description:	Dengue 3- 22kDa is a dengue antigen specially designed for ELISA test, this recombinant peptide contains important epitopes for dengue IgG and IgM antibody recognition. Its applications have been published by the journal. This protein is fused with 6-His fusion partner.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>95% pure (10% SDS-PAGE, Coomassie blue stain). Proprietary chromatographic technique
Concentration:	1.06mg/ml (lot 04192024)
Buffer:	Phosphate buffered saline, pH 7.4
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.



Dengue 3 envelope
22kDa (N-terminal),
10% SDS-PAGE test
on 04-18-2024