

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	Den4-L
Description:	Large size recombinant Dengue 4 envelope antigen is a 45 kDa peptide from 43 to 413 amino acids, this recombinant peptide contains important epitopes for dengue antibody recognition. This protein is fused with 6-His fusion partner at its C-terminus.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>95% pure 12% PAGE (Coomassie). Proprietary chromatographic technique
Concentration:	1.3mg/ml (lot 11-17-2022)
Buffer:	PBS with 25mM K ₂ CO ₃
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

