

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K

Certificate of Analysis

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	D4- NS1
Description:	Recombinant Dengue 4 NS1 is a full length NS1 antigen, it contains multiple epitopes for dengue antibody recognition and vaccine development, it migrated at 45 kDa on 12% SDS-PAGE gel, this protein is fused with 6-His fusion partner.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>95% pure (12% SDS-PAGE, Coomassie blue stain). Proprietary chromatographic technique
Concentration:	0.75mg/ml (Lot 10-13-2013)
Buffer:	Phosphate buffered saline and 50mM arginine
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

