


Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Product Specification Sheet

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	D2-preM
Description:	Dengue virus subtype 2 full length pre-membrane is expressed from <i>E. coli</i> , this protein has a molecular weight 18kDa, it is considered as the target for vaccine development and potential antigen for the diagnosis of dengue fever infection. Recombinant dengue 2 pre-membrane is fused with 6-His tag at its C-terminal.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>95% pure (10% PAGE-PAGE, Coomassie blue stain). Proprietary chromatographic technique
Concentration:	1.3mg/ml (101624)
Buffer:	PBS and 0.25mM K ₂ CO ₃
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.



SDS-PAGE
10-16-2024