

Bioassayplus Inc.

4510 Rhodes Dr. Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Product name:	Recombinant hepatitis B virus core antigen
Catalog #:	rHBc
Description:	Recombinant hepatitis B virus core antigen is a 23kDa peptide showed on 12% SDS-PAGE, core antigen sequence is derived from HBV <i>adv2</i> core region from 1 to185 amino acids, this peptide is fused with 6 x His fusion partner.
Source:	<i>E. coli</i>
Format:	Purified, Liquid
Purification:	>90% pure (10% SDS-PAGE, Coomassie blue stain).
Concentration:	1.08 mg/ml (lot: 11/2019)
Buffer:	25mM Tris base, 10mM K ₂ CO ₃
Applications:	Immunoassay
Storage:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

