

# Bioassayplus Inc.

4510 Rhodes Dr. Windsor, Canada N8W 5k5

## Certificate of Analysis

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Product name:** Recombinant hepatitis C virus core antigen

**Catalog number:** HCV2b-core

**Description:** Recombinant hepatitis C virus genotype 2b core protein contains 120 amino acids, a GST tag is linked to this protein at its N terminus. Recombinant HCV core has a molecular weight 38 kDa.

**Source:** *E. coli*

**Format:** Purified, Liquid

**Purification:** >95% pure (10% SDS-PAGE Coomassie blue staining).

**Concentration** 2.1mg/ml ( lot: E04-20-2022 )

**Buffer:** 6M urea in PBS/25mM K<sub>2</sub>CO<sub>3</sub>

**Applications:** Immunoassay

**Storage:** Upon receipt, aliquot into 1.5ml vials, store at -20°C



GST-HCV2b core, 04-20-2022