

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog # HPV18-LC

Product name: Recombinant human papillomavirus 18 large capsid

Description: Recombinant rHPV18 large capsid is a 55 kDa protein covering full length of HPV18 large capsid, its N-terminus is fused with GST tag. The total molecular size is 78 kDa. HPV18 large capsid is usually used for vaccine study.

Source: *E. coli*

Format: Purified, Liquid

Purification: >90% pure on 10% SDS-PAGE, Coomassie blue stain). Proprietary chromatographic technique

Concentration: 0.72mg/ml (lot 06-4-2018)

Buffer: PBS with 2M urea

Applications: Immunoassay

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

