

Bioassayplus Inc.

4510 Rhodes Drive Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: HPV6-LC

Product name: Recombinant human papillomavirus 6 large capsid

Description: Recombinant rHPV6 antigen is a 55.60 kDa protein covering full length of HPV6 large capsid, its N-terminus is fused with GST tag, the total expressed molecular size is about 78kDa.

Source: *E. coli*

Format: Purified, Liquid

Purification: >90% pure on 10% SDS-PAGE (Coomassie blue stain). Proprietary chromatographic technique

Concentration: 0.9mg/ml (03-2014, requantitated concentration)

Buffer: PBS containing 3M urea and 50mM arginine

Applications: Immunoassay

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.



