

# Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

## Product Specification sheet

<b>Important Note</b>	Centrifuge before opening to ensure complete recovery of vial contents.
<b>Catalog #:</b>	poly-lassa NP
<b>Description:</b>	Recombinant full length Lassa virus nucleoprotein is immunized to BALB/c mouse, mouse IgG is extracted from the mouse serum by protein A column, the purified mouse IgG to Lassa NP protein strongly reacts to recombinant Lassa virus nucleoprotein by ELISA assay.
<b>Source:</b>	Mouse
<b>Format:</b>	90% purified, Liquid
<b>Purification:</b>	Purified by protein A chromatography.
<b>Concentration:</b>	1.5mg/ml
<b>Buffer:</b>	10mM Phosphate buffered saline pH 7.2
<b>Preservative:</b>	0.025% sodium azide
<b>Applications:</b>	Immunoassay
<b>Storage:</b>	Upon receipt, store at 4°C, or -20°C for long time storage