

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, ON, Canada N8W 5K5

Product specification sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog : Mab-mumps

Description: Monoclonal antibody to Mumps virus antibody is specific for the nucleoprotein of Mumps virus (see below). Functions in IFA & ELISA.

Source: mouse

Format: purified IgG

Subtype 1: IgG2a

Purification: purified by protein A chromatography with over 95% purity

Concentration: 3.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: IFA, ELISA

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing