

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Product Specification Sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog: Mab-CP

Description: Specific for *C. pneumoniae* (TWAR)*. Functions in ELISA & IFA. Not reactive with *Chlamydia trachomatis*, *Klebsiella pneumoniae*, Influenza A, Influenza B, Adenovirus, *Mycoplasma pneumoniae*, Parainfluenza type 1, Parainfluenza type 2, Parainfluenza type 3, Respiratory Syncytial virus, *Legionella pneumoniae*, *Haemophilus influenzae* type B, *Bordetella pertussis*, *Streptococcus pneumoniae*, and group A *Streptococcus*

Source: Mouse

Format: Purified IgG

Subtype 1: IgG2b

Purification: Purified by protein A chromatography with over 95% purity

Concentration: 3.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: Immunoassay

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing