

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, ON, Canada N8W 5K5

Product specification sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog: Mab-Trich

Description: Specific to the p65 adhesive antigen of *T. vaginalis*. Functions in immunofluorescence (IFA) & ELISA

Source: Mouse

Format: Purified IgG

Subtype: IgG1

Purification: Purified by protein A chromatography with over 95% purity

Concentration: 3.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: Immunoassay, lateral flow immunoassay

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing