

# Bioassayplus Inc.

4510 Rhodes Dr., Windsor, ON, Canada N8W 5K5

## Product Specification Sheet

<b>Important Note:</b>	Centrifuge before opening to ensure complete recovery of vial contents.
<b>Catalog #:</b>	Pab-den1-NS1
<b>Description:</b>	Polyclonal antibody to dengue-1 NS1 IgG was purified from rabbit serum immunized with whole length recombinant dengue NS1, reactive to all for dengue serotype NS1.
<b>Source:</b>	Rabbit
<b>Format:</b>	Purified IgG
<b>Purification:</b>	>95% pure (12% SDS-PAGE and Coomassie staining), purified by protein A chromatography.
<b>Concentration:</b>	3.2 mg/ml
<b>Buffer:</b>	0.2M glycine pH7.5
<b>Preservative:</b>	0.02% NaN <sub>3</sub>
<b>Applications:</b>	Immunoassay
<b>Storage:</b>	Upon receipt, store at -4°C, and -20°C for long time storage , avoid multiple freeze/thaw cycles.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**