

BioassayPlus Inc.

4510 Rhodes Dr. Windsor, Canada N8W 5K5

Product Specification Sheet

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	Pab-den3-NS1
Description:	Polyclonal antibody to dengue-3 NS1 IgG was derived from the rabbit immunized with whole length recombinant dengue NS1, reactive to all other dengue serotype NS1.
Source:	Rabbit
Format:	Purified IgG
Purification:	>95% pure (12% SDS-PAGE and Coomassie staining), purified by protein A chromatography.
Concentration:	2.9 mg/ml
Buffer:	0.2M glycine, pH7.5
Preservative:	0.02% sodium azide
Applications:	Immunoassay
Storage:	Upon receipt, store at -4°C, and -20°C for long time storage , avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.