

# Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

## Product Specification Sheet

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Product Name :** Goat antibody to Chlamydia trachomatis

**Catalog:** Pab-C. trachomatis

**Description:** Goat IMMUNOGEN : Purified MOMP from strain L2, SPECIFICITY : Major Outer Membrane Protein (MOMP) TITER : >1:1,000 by indirect immunofluorescence vs. all serovars (A-K, L1-L3) CROSS REACTIVITY: Does not react with C. psittacii or C. pneumoniae in MIF UNINFECTED CELL REACTIVITY: Negative vs. HEp-2 cells and egg yolk sac

**Source:** Goat

**Format:** Purified IgG

**Purification:** Purified by protein A chromatography with over 95% purity

**Concentration:** 4.0mg/ml

**Buffer:** PBS

**Preservative:** 0.05% Sodium nitrate

**Applications:** Immunoassay

**Storage:** Upon receipt, store at -4°C. Avoid multiple freezing/thawing