

# Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

## Product specification sheet

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog :** Mab-SLT-1B

**Description:** Mab-SLT-1B is specific for Shiga toxin type 1- B subunit (Stx-1b). Functions well in ELISA

**Source:** mouse

**Format:** purified IgG

**Subtype 1:** IgG1

**Purification:** purified by protein A chromatography with over 95% purity

**Concentration:** 3.0mg/ml

**Buffer:** PBS

**Preservative:** 0.05% Sodium nitrate

**Applications:** Lateral flow immunoassay

**Storage:** Upon receipt, store at -4°C. Avoid multiple freezing/thawing