

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, Canada N8W 5K5

Product specification sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog : Mab-SLT-2B

Description: Mab-SLT-2B is specific for Shiga toxin type 2 - B subunit (Stx-2b). Functions well in ELISA

Source: mouse

Format: Purified IgG

Subtype 1: IgG1

Purification: purified by protein A chromatography with over 95% purity

Concentration: 2.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: Immunoassay

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing