

Bioassayplus Inc.

4510 Rhodes Dr. Windsor, Canada N8W 5K5

Certificate of Analysis

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product name: Recombinant Ureaplasma Urealyticum antigen (strain serovar 3)

Catalog # : UU-3MB

Description: Ureaplasma Urealyticum Serovar 3 is a common pathogen to cause sexual transmitted disease. MB antigen peptide fused with GST from Ureaplasma Urealyticum Serovar 3 is expressed from *E. coli*, the expressed fusion protein migrates at 40kDa position on SDS-PAGE.

Source: *E. coli*

Format: Purified, Liquid


Purification: >90% pure (10% SDS-PAGE, Coomassie blue stain). Proprietary chromatographic technique

Concentration: 4.58mg/ml (lot 05-2017)

Buffer: Phosphate buffered saline and 50mM arginine

Applications: Immunoassay

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.



April 20, 2022, SDS-PAGE
test for Recombinant
Ureaplasma Urealyticum
antigen (strain serovar 3)