

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, ON, Canada N8W 5K5

Product specification sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog : Mab-WNV-NS1

Description: WNV belongs to the genus of Flavivirus family. WNV is a positive sense RNA virus, encoding seven nonstructural proteins and three structural proteins. West Nile fever infection spreads by mosquitoes.

West Nile virus monoclonal antibody is specific for WNV envelope for ELISA and other immunoassay purpose.

Source: mouse

Format: purified IgG

Subtype 1: IgG1

Purification: purified by protein A chromatography with over 95% purity

Concentration: 3.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: Lateral flow immunoassay

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing