

**Antigen conjugated gold supplied either by 10 gold pads or 5ml 5OD gold solution**

<b>Gold pad made for disease</b>	<b>Catalog</b>	<b>Conjugated antigen</b>	<b>Application</b>
ASFV (Africa swine fever virus)	G-ASFV	ASFV-P72	Antibody test
CHFV (Congo hemorrhagic fever)	G-CHFV	CHFV-NP	Antibody test
Covid 19	G-covid	Covid19 NP	Antibody test
Ebola Z virus	G-EBO	Ebo-Z VP40	Antibody test
Dengue virus	G-Den	Den epi-10 / 13	Antibody test
Hepatitis A virus	G-HAV	HAV-VP3	Antibody test
Hepatitis B	G-HBe	HBe	HBe test
Hepatitis C virus	G-HCV	HCV cocktail	Antibody test
Human immunodeficiency virus	G-HIV	HIV1 gp41 / HIV2 gp36	Antibody test
Japanese encephalitis virus	G-JEV	JEV envelope	Antibody test
Rubella virus	G-Rubella	Rubellar capsid	Antibody test
Zika virus	G-zika	Zika envelope	Antibody test
Scrub Typhus	G-scrube typhus	O. tsutsugamushi r56	Antibody test
Toxoplasma Gondii	G-SAG2	Toxoplasma SAG2	Antibody test
	G-GRA6	Toxoplasma GRA6	Antibody test
	G-GRA7	Toxoplasma GRA7	Antibody test
	G-GRA8	Toxoplasma GRA8	Antibody test
T. Pallidum ( Syphilis )	G1-syphilis	T. Pallidum Tp17 / TP47	Antibody test
	G2-syphilis	Mosaic TP and Tp17	Antibody test
Flariasis	G-Flariasis	Wb123	Antibody test
Leshimania	G-leish	Leishmania k39	Antibody test
Malaria	G-malaria	Pf and Pv MSP1	Antibody test
TSH receptor	G-TSH recptor	TSH receptor	Antibody to TSH receptor