

Bioassayplus Inc.

4510 Rhodes Dr. Windsor, ON, Canada N8W 5K5

Product Specification Sheet

- Product name:** Gold conjugated *Streptavidin* (*Streptomyces avidinii*)
- Catalog #:** Strep-gold
- Description:** Gold conjugated *Streptavidin* is nature streptavidin conjugated to 40nm colloid gold, which can be applied to gold pad with the target conjugate together for the control line in the development of lateral flow assay. Biotin-BSA is available in our company used to coat to the cellulose membrane making control line.
- Concentration:** OD5.0
- Expiration date:** One year after receiving the product
- Storage :** 4°C for long term
- Appearance:** Show homogeneous red solution without precipitation, sometimes if show small particle, Use poppet up and down gently then use it.