

Bioassayplus Inc.

4510 Rhodes Dr., Windsor, ON, Canada N8W 5K5

Product specification sheet

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog : Mab-M. genitalium

Description: Specific for M. genitalium. Specific to M. genitalium. Functions in ELISA (pairing below) & IFA.
Non-reactive with: Mycoplasma hominis, Mycoplasma pneumoniae, Ureaplasma urealyticum,
Herpes virus type 2, Chlamydia trachomatis, Neisseria gonorrhoeae, and human cells

Source: Mouse

Format: Purified IgG

Subtype 1: IgG2b

Purification: purified by protein A chromatography with over 95% purity

Concentration: 3.0mg/ml

Buffer: PBS

Preservative: 0.05% Sodium nitrate

Applications: Immunoassay

Storage: Upon receipt, store at -4°C. Avoid multiple freezing/thawing