

Certificate of Analysis

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product name: Mosaic recombinant cytomegalovirus antigen

Catalog # : rCMV

Description: Recombinant cytomegalovirus antigens contains the multiple epitopes from P150-P52-P38-p65-P28, the molecular weight is about 40kDa, this recombinant protein is fused with 6 x his tag.

Source: *E. coli*

Format: Purified, Liquid

Purification: >95% pure (12% PAGE Coomassie blue stain). Proprietary chromatographic technique

Concentration: 0.90mg/ml (Lot: 10/18/2022)

Buffer: Phosphate buffered saline/25mM K₂CO₃ with 2M urea

Applications: Immunoassay

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.



New SDS-PAGE
test, 09-21-2024